

REG4 antibody

Catalog Number: 114683

Product name

REG4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

REG4 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in SGC-7901 cells, mouse liver tissue. Observed molecular weight by Western-blot: 18 kDa

Preparation

This antibody was obtained by immunization of REG4 recombinant protein (Accession Number: NM_032044). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IHC: 1:50-1:500

Validations



SGC-7901 cells were subjected to SDS PAGE followed by western blot with Catalog No:114683(REG4 antibody) at dilution of 1:300 

Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:114683(REG4 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:114683(REG4 Antibody) at dilution of 1:200 (under 40x lens).