

RIC8A antibody



Catalog Number: 114698

Product name

RIC8A antibody

Specificity

Human; other species not tested.

Antibody description

RIC8A Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Positive IHC detected in human brain tissue, human pancreas cancer tissue. Observed molecular weight by Western-blot: 60-70 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: BC011821). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114698(RIC8A antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:114698(RIC8A Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:114698(RIC8A Antibody) at dilution of 1:50 (under 40x lens)