

# RILPL1 antibody



Catalog Number: 114701

## Product name

RILPL1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RILPL1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of RILPL1 recombinant protein (Accession Number: NM\_001319244). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114701(RILPL1 Antibody) at dilution of 1:600