

# RIP3 antibody

Catalog Number: 114714

## Product name

RIP3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RIP3 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, NIH/3T3 cells. Positive IP detected in SW 1990 cells. Observed molecular weight by Western-blot: 57kd

## Preparation

This antibody was obtained by immunization of RIP3 recombinant protein (Accession Number: NM\_006871). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:114714(RIPK3 Antibody) at dilution of 1:300



IP Result of anti-RIPK3 (IP:Catalog No:114714, 3ug; Detection:Catalog No:114714 1:300) with SW 1990 cells lysate 3000ug.