

RNF141 antibody

Catalog Number: 114741

Product name

RNF141 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

RNF141 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 26 kDa

Preparation

This antibody was obtained by immunization of RNF141 recombinant protein (Accession Number: NM_016422). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114741(RNF141 Antibody) at dilution of 1:300



IP Result of anti-RNF141 (IP:Catalog No:114741, 4ug; Detection:Catalog No:114741 1:300) with mouse testis tissue lysate 4000ug.