

# RNGTT antibody



Catalog Number: 114766

## Product name

RNGTT antibody

## Immunogen

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RNGTT Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western-blot: 69 kDa

## Preparation

This antibody was obtained by immunization of RNGTT recombinant protein (Accession Number: NM\_003800). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:114766(RNGTT antibody) at dilution of 1:800



IP Result of anti-RNGTT (IP:Catalog No:114766, 4ug; Detection:Catalog No:114766 1:500) with mouse kidney tissue lysate 4000ug.