

# RNPC3 antibody

Catalog Number: 114770

## Product name

RNPC3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RNPC3 Rabbit Polyclonal antibody. Positive IP detected in mouse spleen tissue. Positive WB detected in mouse spleen tissue, 3T3-L1 cells. Observed molecular weight by Western-blot: 62-65 kDa

## Preparation

This antibody was obtained by immunization of RNPC3 recombinant protein (Accession Number: NM\_017619). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:114770(RNPC3 Antibody) at dilution of 1:300



IP Result of anti-RNPC3 (IP:Catalog No:114770, 4ug; Detection:Catalog No:114770 1:300) with mouse spleen tissue lysate 4000ug.