

# ROBO3 antibody

Catalog Number: 114774

## Product name

ROBO3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ROBO3 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in mouse stomach tissue. Observed molecular weight by Western-blot: 70kd

## Preparation

This antibody was obtained by immunization of ROBO3 recombinant protein (Accession Number: BC008623). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



mouse stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:114774(ROBO3 antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:114774(ROBO3 antibody) at dilution of 1:50 (under 10x lens)