

# ROBO3-Specific antibody



Catalog Number: 114775

## Product name

ROBO3-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ROBO3-Specific Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue, human brain tissue. Positive WB detected in mouse brain tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 55kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_022370). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114775(ROBO3-Specific antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human gliomas slide using Catalog No:114775(ROBO3-Specific Antibody) at dilution



of 1:50

Immunohistochemistry of paraffin-embedded human gliomas slide using Catalog No:114775(ROBO3-Specific Antibody) at dilution of 1:50