

RORA antibody



Catalog Number: 114788

Product name

RORA antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RORA Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HeLa cells, Jurkat cells, K-562 cells, mouse kidney tissue, PC-3 cells, rat liver tissue. Observed molecular weight by Western-blot: 70-75kd

Preparation

This antibody was obtained by immunization of RORA recombinant protein (Accession Number: NM_134261). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114788(RORA antibody) at dilution of 1:1000