

RXFP4 antibody



Catalog Number: 114944

Product name

RXFP4 antibody

Specificity

Human; other species not tested.

Antibody description

RXFP4 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_181885).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:114944(RXFP4 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:114944(RXFP4 Antibody) at dilution of 1:200 (under 40x lens).