

SCLT1 antibody

Catalog Number: 114999

Product name

SCLT1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SCLT1 Rabbit Polyclonal antibody. Positive IP detected in human placenta tissue. Positive WB detected in mouse cerebellum tissue. Observed molecular weight by Western-blot: 75-80 kDa

Preparation

This antibody was obtained by immunization of SCLT1 recombinant protein (Accession Number: NM_001300897). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with Catalog No:114999(SCLT1 antibody) at dilution of 1:300



IP Result of anti-SCLT1 (IP:Catalog No:114999, 4ug; Detection:Catalog No:114999 1:1000) with human placenta tissue lysate 4000ug.