

SAMSN1 antibody



Catalog Number: 115058

Product name

SAMSN1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SAMSN1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 42-50 kDa

Preparation

This antibody was obtained by immunization of SAMSN1 recombinant protein (Accession Number: NM_022136). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:115058(SAMSN1 Antibody) at dilution of 1:1500



Immunofluorescent analysis of HepG2 cells using Catalog No:115058(SAMSN1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)