

SEC22B antibody



Catalog Number: 115077

Product name

SEC22B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SEC22B Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells, rat liver tissue. Observed molecular weight by Western-blot: 24 kDa

Preparation

This antibody was obtained by immunization of SEC22B recombinant protein (Accession Number: NM_004892). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:115077(SEC22B antibody) at dilution of 1:500