

Secretogranin II/Chromogranin C antibody



Catalog Number: 115095

Product name

Secretogranin II/Chromogranin C antibody

WB: 1:500-1:5000

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

Secretogranin II/Chromogranin C Rabbit Polyclonal antibody. Positive IHC detected in human cerebellum tissue, human brain tissue, human pancreas tissue, human small intestine tissue. Positive WB detected in SH-SY5Y cells, Y79 cells. Observed molecular weight by Western-blot: 66 kDa



Preparation

This antibody was obtained by immunization of Secretogranin II/Chromogranin C recombinant protein (Accession Number: NM_003469). Purification method: Antigen affinity purified.

SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:115095(SCG2 antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemistry of paraffin-embedded human cerebellum tissue slide using Catalog No:115095(SCG2 Antibody) at dilution of 1:50 (under 10x lens)

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:



Immunohistochemistry of paraffin-embedded

Secretogranin II/Chromogranin C antibody



Catalog Number: 115095

human cerebellum tissue slide using Catalog

No:115095(SCG2 Antibody) at dilution of 1:50
(under 40x lens)