

SEPT4 antibody

Catalog Number: 115118

Product name

SEPT4 antibody

Specificity

Human; other species not tested.

Antibody description

SEPT4 Rabbit Polyclonal antibody. Positive WB detected in A375 cells. Positive IP detected in A375 cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of SEPT4 recombinant protein (Accession Number: NM_080416). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:115118(SEPT4 antibody) at dilution of 1:300



IP Result of anti-SEPT4 (IP:Catalog No:115118, 4ug; Detection:Catalog No:115118 1:300) with A375 cells lysate 2200ug.