

SFRS8 antibody

Catalog Number: 115153

Product name

SFRS8 antibody

Specificity

Human; other species not tested.

Antibody description

SFRS8 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 130-150 kDa

Preparation

This antibody was obtained by immunization of SFRS8 recombinant protein (Accession Number: NM_004592). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115153(SFRS8 Antibody) at dilution of 1:2000



Immunofluorescent analysis of HeLa cells using Catalog No:115153(SFRS8 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)