

SGEF antibody



Catalog Number: 115164

Product name

SGEF antibody

Specificity

human; other species not tested.

Antibody description

SGEF Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells. Observed molecular weight by Western-blot: 97kd

Preparation

This antibody was obtained by immunization of SGEF recombinant protein (Accession Number: BC016628). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115164(SGEF antibody) at dilution of 1:300