

# SH3BP5 antibody



Catalog Number: 115190

## Product name

SH3BP5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SH3BP5 Rabbit Polyclonal antibody. Positive WB detected in Raji cells. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 47kd

## Preparation

This antibody was obtained by immunization of SH3BP5 recombinant protein (Accession Number: NM\_004844). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:115190(SH3BP5 antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:115190(SH3BP5 antibody) at dilution of 1:50 (under 10x lens)