

SIRP alpha antibody



Catalog Number: 115238

Product name

SIRP alpha antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SIRP alpha Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, human heart tissue, mouse cerebellum tissue, PC-12 cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of SIRP alpha recombinant protein (Accession Number: NM_001330728). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:115238(SIRPA Antibody) at dilution of 1:600