

SIRP gamma antibody



Catalog Number: 115240

Product name

SIRP gamma antibody

Specificity

Human; other species not tested.

Antibody description

SIRP gamma Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 55kd

Preparation

This antibody was obtained by immunization of SIRP gamma recombinant protein (Accession Number: NM_080816). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115240(SIRPG antibody) at dilution of 1:400