

SLAMF7 antibody



Catalog Number: 115296

Product name

Polyclonal

SLAMF7 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

SLAMF7 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 37 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of SLAMF7 recombinant protein (Accession Number: NM_001282592). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115296(SLAMF7 antibody) at dilution of 1:300