

# SLC22A18AS antibody



Catalog Number: 115315

## Product name

SLC22A18AS antibody

## Specificity

Human; other species not tested.

## Antibody description

SLC22A18AS Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Observed molecular weight by Western-blot: 30 kDa;65 kDa

## Preparation

This antibody was obtained by immunization of SLC22A18AS recombinant protein (Accession Number: XM\_017017836). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:1000

## Validations



IP Result of anti-SLC22A18AS (IP:Catalog No:115315, 3ug; Detection:Catalog No:115315 1:300) with HepG2 cells lysate 3600ug.