

SLITRK3 antibody

Catalog Number: 115358

Product name

SLITRK3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SLITRK3 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 120-125 kDa

Preparation

This antibody was obtained by immunization of SLITRK3 recombinant protein (Accession Number: NM_001318811). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115358(SLITRK3 antibody) at dilution of 1:200



IP Result of anti-SLITRK3 (IP:Catalog No:115358, 4ug; Detection:Catalog No:115358 1:300) with mouse brain tissue lysate 4000ug.