

SMN antibody



Catalog Number: 115394

Product name

SMN antibody

Specificity

Human; other species not tested.

Antibody description

SMN Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in K-562 cells, HEK-293 cells, HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 38 kDa

Preparation

This antibody was obtained by immunization of SMN recombinant protein (Accession Number: NM_022875). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:115394(SMN2 antibody) at dilution of 1:2000



Immunofluorescent analysis of HepG2 cells using Catalog No:115394(SMN2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)