

SLU7 antibody

Catalog Number: 115413

Product name

SLU7 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

SLU7 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, K-562 cells, Y79 cells, mouse colon tissue. Positive IHC detected in human colon cancer. Observed molecular weight by Western-blot: 68kd

Preparation

This antibody was obtained by immunization of SLU7 recombinant protein (Accession Number: XM_011534418). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115413(SLU7 antibody) at dilution of 1:100



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:115413(SLU7 antibody) at dilution of 1:50 (under 10x lens)