

SMARCD2 antibody



Catalog Number: 115434

Product name

SMARCD2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

SMARCD2 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells, mouse thymus tissue. Positive IHC detected in human cervical cancer. Observed molecular weight by Western-blot: 52kd

Preparation

This antibody was obtained by immunization of SMARCD2 recombinant protein (Accession Number: BC018953). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



RAW264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:115434(SMARCD2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:115434(SMARCD2 antibody) at dilution of 1:50 (under 10x lens)