

SOX21 antibody



Catalog Number: 115516

Product name

SOX21 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SOX21 Rabbit Polyclonal antibody. Positive WB detected in mouse embryo tissue, Jurkat cells. Observed molecular weight by Western-blot: 25-30 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_007084). Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse embryo tissue were subjected to SDS PAGE followed by western blot with Catalog No:115516(SOX21 Antibody) at dilution of 1:300