

# SOX8 antibody



Catalog Number: 115520

## Product name

SOX8 antibody

## Specificity

Human; other species not tested.

## Antibody description

SOX8 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, COLO 320 cells. Positive IHC detected in human brain. Observed molecular weight by Western-blot: 47kd

## Preparation

This antibody was obtained by immunization of SOX8 recombinant protein (Accession Number: NM\_014587). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115520(SOX8 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human brain using Catalog No:115520(SOX8 antibody) at dilution of 1:50 (under 10x lens)