

SP2 antibody

Catalog Number: 115525

Product name

SP2 antibody

Specificity

Human; other species not tested.

Antibody description

SP2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, K-562 cells. Observed molecular weight by Western-blot: 72 kDa

Preparation

This antibody was obtained by immunization of SP2 recombinant protein (Accession Number: XM_047436577). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115525(SP2



Antibody) at dilution of 1:300

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:115525(SP2 Antibody) at dilution of 1:300