

SPANXA2/SPANXA1 antibody



Catalog Number: 115539

Product name

SPANXA2/SPANXA1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human; other species not tested.

Validations

Antibody description

SPANXA2/SPANXA1 Rabbit Polyclonal antibody. Positive IHC detected in human skin cancer tissue, human testis tissue. Positive IF detected in HEK-293 cells. Positive WB detected in DU 145 cells, HEK-293 cells, PC-3 cells. Observed molecular weight by Western-blot: 27 kDa



Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_145662). Purification method: Antigen affinity purified.

DU 145 cells were subjected to SDS PAGE followed by western blot with Catalog No:115539(SPANXA2 antibody) at dilution of

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IF, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000



1:500

Immunofluorescent analysis of HEK-293 cells using Catalog No:115539(SPANXA2 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG




Immunohistochemistry of paraffin-embedded

SPANXA2/SPANXA1 antibody



Catalog Number: 115539

human testis tissue slide using Catalog
No:115539(SPANXA2 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded
human skin cancer tissue slide using Catalog
No:115539(SPANXA2 Antibody) at dilution of 1:50