

# SPO11 antibody



Catalog Number: 115576

## Product name

SPO11 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

SPO11 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human testis tissue. Observed molecular weight by Western-blot: 55-70 kDa

## Preparation

This antibody was obtained by immunization of SPO11 recombinant protein (Accession Number: NM\_198265). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115576(SPO11 antibody) at dilution of 1:400