

SSU72 antibody



Catalog Number: 115616

Product name

SSU72 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SSU72 Rabbit Polyclonal antibody. Positive WB detected in 293/HeLa cells.

Preparation

This antibody was obtained by immunization of SSU72 recombinant protein (Accession Number: NM_014188). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: N/A

Validations



WB result of SSU72 antibody (Catalog No:115616, 1:1,000) with 293, HeLa cell lysate by Dr. Katherine Jones' lab in the Salk Institute for Biological Studies.