

# ST7L antibody



Catalog Number: 115632

## Product name

ST7L antibody

## Immunogen

## Specificity

Human,Mouse,Rat; other species not tested.

## Antibody description

ST7L Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 60-66 kDa

## Preparation

This antibody was obtained by immunization of ST7L recombinant protein (Accession Number: NM\_017744). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115632(ST7L antibody) at dilution of 1:1500



IP Result of anti-ST7L (IP:Catalog No:115632, 3ug; Detection:Catalog No:115632 1:1500) with mouse testis tissue lysate 5000ug.