

# Sprr2i antibody



Catalog Number: 115646

## Product name

Sprr2i antibody

## Specificity

Human; other species not tested.

## Antibody description

Sprr2i Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue.

## Preparation

This antibody was obtained by immunization of Sprr2i recombinant protein (Accession Number: NM\_011475). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions

Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:115646(Sprr2i Antibody) at dilution of 1:50 (under 40x lens)