

SPZ1 antibody

Catalog Number: 115657

Product name

SPZ1 antibody

Specificity

human; other species not tested.

Antibody description

SPZ1 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, human testis tissue, DU 145 cells. Positive IHC detected in human testis. Observed molecular weight by Western-blot: 49-55kd

Preparation

This antibody was obtained by immunization of SPZ1 recombinant protein (Accession Number: NM_032567). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115657(SPZ1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human testis using Catalog No:115657(SPZ1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:115657(SPZ1 antibody) at dilution of 1:100 (under 40x lens)