

STRA6 antibody



Catalog Number: 115733

Product name

STRA6 antibody

Specificity

Human; other species not tested.

Antibody description

STRA6 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 95 kDa

Preparation

This antibody was obtained by immunization of STRA6 recombinant protein (Accession Number: NM_022369). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115733(STRA6 Antibody) at dilution of 1:300