

TARS antibody

Catalog Number: 115853

Product name

TARS antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TARS Rabbit Polyclonal antibody. Positive WB detected in A549 cells. Positive IF detected in MCF-7 cells, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 66-72kd

Preparation

This antibody was obtained by immunization of TARS recombinant protein (Accession Number: NM_152295). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:115853(TARS antibody) at dilution of 1:600



Immunofluorescent analysis of MCF-7 cells, using TARS antibody Catalog No:115853 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).