

TBC1D24 antibody



Catalog Number: 115867

Product name

TBC1D24 antibody

Immunogen

[Human TBC1D24 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

TBC1D24 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 60-63 kDa

Preparation

This antibody was obtained by immunization of TBC1D24 recombinant protein (Accession Number: XM_017023493). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115867(TBC1D24 Antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells using Catalog No:115867(TBC1D24 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)