

TCEB1 antibody

Catalog Number: 115900

Product name

TCEB1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TCEB1 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in mouse testis tissue, HeLa cells. Observed molecular weight by Western-blot: 12 kDa

Preparation

This antibody was obtained by immunization of TCEB1 recombinant protein (Accession Number: XM_047422155). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115900(TCEB1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115900(TCEB1 antibody) at dilution of 1:50 (under 25x lens)