

TCP10 antibody

Catalog Number: 115916

Product name

TCP10 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TCP10 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, HepG2 cells, mouse testis tissue. Observed molecular weight by Western-blot: 46kd

Preparation

This antibody was obtained by immunization of TCP10 recombinant protein (Accession Number: NM_001366317). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115916(TCP10 antibody) at dilution of 1:500