

TEAD2 antibody



Catalog Number: 115932

Product name

Polyclonal

TEAD2 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

TEAD2 Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue, mouse lung tissue. Observed molecular weight by Western-blot: 55-65 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of TEAD2 recombinant protein (Accession Number: XM_047439527). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:115932(TEAD2 Antibody) at dilution of 1:300

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality