

TEAD4 antibody

Catalog Number: 115934

Product name

TEAD4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TEAD4 Rabbit Polyclonal antibody. Positive IP detected in MCF-7 cells. Positive WB detected in MDA-MB-453s cells, BGC-823 cells, MCF-7 cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of TEAD4 recombinant protein (Accession Number: NM_003213). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:115934 (TEAD4 Antibody) at dilution of 1:600



IP Result of anti-TEAD4 (IP: Catalog No:115934, 3ug; Detection: Catalog No:115934 1:300) with MCF-7 cells lysate 2000ug.