

TFAP2A,AP-2 antibody



Catalog Number: 115998

Product name

TFAP2A,AP-2 antibody

Immunogen

[Human TFAP2A,AP-2 Recombinant protein \(His tag\)](#)

Specificity

Human,Mouse,Rat; other species not tested.

Antibody description

TFAP2A,AP-2 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of TFAP2A,AP-2 recombinant protein (Accession Number: NM_001032280). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog

No:115998(TFAP2A,AP-2 antibody) at dilution of 1:300