

# TFAP2B antibody

Catalog Number: 115999

## Product name

TFAP2B antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TFAP2B Rabbit Polyclonal antibody. Positive IP detected in MCF-7 cells. Positive WB detected in MCF7 cells. Observed molecular weight by Western-blot: 50kd

## Preparation

This antibody was obtained by immunization of TFAP2B recombinant protein (Accession Number: NM\_003221). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:115999(TFAP2B antibody) at dilution of 1:800



IP Result of anti-TFAP2B (IP:Catalog No:115999, 4ug; Detection:Catalog No:115999 1:500) with MCF-7 cells lysate 800ug.