

TGIF2LX antibody



Catalog Number: 116027

Product name

TGIF2LX antibody

Specificity

human; other species not tested.

Antibody description

TGIF2LX Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 34kd

Preparation

This antibody was obtained by immunization of TGIF2LX recombinant protein (Accession Number: BC029920). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116027(TGIF2LX antibody) at dilution of 1:600