

# THG1L antibody



Catalog Number: 116043

## Product name

THG1L antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

THG1L Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human endometrial cancer tissue.

## Preparation

This antibody was obtained by immunization of THG1L recombinant protein (Accession Number: BC001852). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human endometrial cancer using Catalog No:116043(THG1L antibody) at dilution of 1:100 (under 10x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116043(THG1L Antibody) at dilution of 1:2000