

# Thrombomodulin antibody



Catalog Number: 116054

## Product name

Thrombomodulin antibody

## Specificity

Human; other species not tested.

## Antibody description

Thrombomodulin Rabbit Polyclonal antibody.  
Positive IHC detected in human heart tissue.  
Positive WB detected in human heart tissue.  
Observed molecular weight by Western-blot: 105 kDa

## Preparation

This antibody was obtained by immunization of Thrombomodulin recombinant protein (Accession Number: NM\_000361). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116054(THBD antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human heart using Catalog No:116054(THBD antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:116054(THBD antibody) at dilution of 1:50 (under 40x lens)