

Thrombomodulin antibody



Catalog Number: 116054

Product name

Thrombomodulin antibody

IHC: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

Thrombomodulin Rabbit Polyclonal antibody.
Positive IHC detected in human heart tissue.
Positive WB detected in human heart tissue.
Observed molecular weight by Western-blot: 105 kDa

Preparation

This antibody was obtained by immunization of Thrombomodulin recombinant protein (Accession Number: NM_000361). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116054(THBD antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human heart using Catalog No:116054(THBD antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:116054(THBD antibody) at dilution of 1:50 (under 40x lens)