

TMCO6 antibody

Catalog Number: 116084

Product name

TMCO6 antibody

WB: 1:1000-1:10000

IHC: 1:20-1:200


Immunogen

[Human TMCO6 Recombinant protein \(GST tag\)](#)

Validations

Specificity

human, mouse, rat; other species not tested.




Antibody description

TMCO6 Rabbit Polyclonal antibody. Positive IHC detected in human spleen. Positive WB detected in human testis tissue, mouse brain tissue. Observed molecular weight by Western-blot: 50-54kd

human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116084(TMCO6 antibody) at dilution of 1:500

Preparation

This antibody was obtained by immunization of TMCO6 recombinant protein (Accession Number: XM_024446125). Purification method: Antigen affinity purified.



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

Immunohistochemical of paraffin-embedded human spleen using Catalog No:116084(TMCO6 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human spleen using Catalog No:116084(TMCO6 antibody) at dilution of 1:100 (under 10x lens)

TMCO6 antibody



Catalog Number: 116084

antibody) at dilution of 1:100 (under 40x lens)